

## Gallios

Laser	Channel	BP	Fluorochromes
Violet (405 nm)	FL10	550/40	<ul style="list-style-type: none"> <li>• Pacific Orange</li> <li>• Krome Orange</li> <li>• BV570</li> <li>• Qdot 565</li> </ul>
	FL9	450/50	<ul style="list-style-type: none"> <li>• Pacific Blue</li> <li>• BV421</li> <li>• efluor 450</li> <li>• CFP</li> <li>• DAPI</li> </ul>
Blue (488 nm)	FL5	755 LP	<ul style="list-style-type: none"> <li>• PE/Cy7</li> </ul>
	FL4	695/30	<ul style="list-style-type: none"> <li>• PE/Cy5</li> <li>• PerCP</li> <li>• 7-AAD</li> <li>• PE/Cy5.5</li> <li>• PerCP/Cy5.5</li> </ul>
	FL3	620/20	<ul style="list-style-type: none"> <li>• PE-TexasRed</li> <li>• PI</li> </ul>
	FL2	575/25	<ul style="list-style-type: none"> <li>• PE</li> <li>• Alexa555</li> </ul>
	FL1	525/40	<ul style="list-style-type: none"> <li>• FITC</li> <li>• GFP</li> <li>• Alexa488</li> <li>• CFSE</li> <li>• MitoTracker Green</li> </ul>
Red (640 nm)	FL8	755 LP	<ul style="list-style-type: none"> <li>• APC/Cy5.5</li> <li>• APC/Cy7</li> <li>• APC/efluor 780</li> <li>• Alexa780</li> <li>• Live-Dead Near IR</li> </ul>
	FL7	725/20	<ul style="list-style-type: none"> <li>• Alexa700</li> </ul>
	FL6	660/20	<ul style="list-style-type: none"> <li>• APC</li> <li>• Alexa647</li> <li>• Cy5</li> <li>• MitoTracker DeepRed</li> </ul>